

GSI Caspase-2 Colorimetric Assay Kit DataSheet

Caspase 2 (CASP2) is an enzyme that, in humans, is encoded by the CASP2 gene.^[1] CASP2 orthologs have been identified in nearly all mammals for which complete genome data are available. Caspase 2 proteolytically cleaves other proteins. It belongs to a family of cysteine proteases called caspases that cleave proteins only at an amino acid following an aspartic acid residue. Within this family, caspase 2 is part of the Ich-1 subfamily. It is one of the most conserved caspases in different species of animal. Caspase 2 has a similar amino acid sequence to initiator caspases, including caspase 1, caspase 4, caspase 5, and caspase 9. It is produced as a zymogen, which contains a long pro-domain that is similar to that of caspase 9 and contains a protein interaction domain known as a CARD domain. Pro-caspase-2 contains two subunits, p19 and p12. It has been shown to associate with several proteins involved in apoptosis using its CARD domain, including RIP-associated Ich-1/Ced-3-homologue protein with a death domain (RAIDD), apoptosis repressor with caspase recruitment domain (ARC), and death effector filament-forming Ced-4-like apoptosis protein (DEFCAP).^[2] Together with RAIDD and p53-induced protein with a death domain ([PIDD])(LRDD), caspase 2 has been shown to form the so-called PIDDosome,^[3] which may serve as an activation platform for the protease, although it may also be activated in the absence of PIDD.^[4] Overall, caspase 2 appears to be a very versatile caspase with multiple functions beyond cell death induction.^{[5][6]}

References

1. Kumar S, et al. (1995). *Hum. Genet.* 95 (6): 641–4.
2. Zhivotovsky B, Orrenius S (2005). *Biochem. Biophys. Res. Commun.* 331 (3): 859–67.
3. Tinel A, Tschopp J (2004). *Science.* 304 (5672): 843–6.
4. Manzl C, et al. (2009). *J. Cell Biol.* 185 (2): 291–303.
5. Krumschnabel G, et al. (2009). *Oncogene.* 28 (35): 3093–6.
6. Krumschnabel G, et al. (2009). *Cell Death Differ.* 16 (2): 195–207.

PRINCIPLE OF THE ASSAY

Caspase 2 has a substrate specificity by recognizing peptide motif x-Asp-x-x-Asp. The C-terminal Asp is absolutely required while variations at other four positions can be tolerated. Caspase-2 cleave a variety of cellular substrates after aspartic acid residues—a characteristic that is central to their role in mammalian apoptosis. The Caspase-2 Colorimetric Assay Kit provides a simple and convenient means for assaying the activity of caspases that recognize the tetrapeptide sequence VDVAD. The assay is based on spectrophotometric detection of the chromophore p-nitroaniline (p-NA) after cleavage from the p-NA-labeled Ac-VDVAD-pNA. The p-NA light emission can be quantified using a spectrophotometer or a microtiter plate reader at 400–410nm. Comparison of the absorbance of p-NA from an apoptotic sample with an uninduced control allows determination of the fold increase in caspase-2 activities.



GSI Caspase-2 Colorimetric Assay Kit DataSheet

This package insert must be read in its entirety before using this product.

Storage

Store at 4°C as supplied for 6 months. Store at -20°C after reconstitution for up to 3 months.

MATERIALS PROVIDED

Description	Quantity	Description	Quantity	Description	Quantity
Substrate	1	2 x Assay Buffer	1	200 mM DTT	1
DataSheet/Manual	1				

Bring all reagents to room temperature before use.

Reagent Preparations

2 x Assay Buffer, 20 mL- Dilute to 1 x Assay Buffer with deionized and distill water prior to use.

Substrate: Reconstitute in 1.1 mL of DMSO and dissolve completely prior to use and store at -20°C.

200 mM DTT: Reconstitute in 0.55 mL of deionized and distill water prior to use and store at -20°C.

For Research Use Only

Related products

GR107005 GSI Caspase-3/7 Colorimetric Assay Kit
GR107006 GSI Caspase-3/7 Fluorometric Assay Kit
GR107007 GSI Caspase-1 Colorimetric Assay Kit
GR107008 GSI Caspase-1 Fluorometric Assay Kit
GR107018 GSI Caspase-2 Fluorometric Assay Kit
GR107019 GSI Caspase-4 Colorimetric Assay Kit
GR107020 GSI Caspase-4 Fluorometric Assay Kit
GR107021 GSI Caspase-5 Colorimetric Assay Kit
GR107022 GSI Caspase-5 Fluorometric Assay Kit
GR107023 GSI Caspase-6 Colorimetric Assay Kit
GR107024 GSI Caspase-6 Fluorometric Assay Kit
GR107025 GSI Caspase-8 Colorimetric Assay Kit
GR107026 GSI Caspase-8 Fluorometric Assay Kit
GR107027 GSI Caspase-9 Colorimetric Assay Kit
GR107028 GSI Caspase-9 Fluorometric Assay Kit
GR107029 GSI Caspase-10 Colorimetric Assay Kit
GR107030 GSI Caspase-10 Fluorometric Assay Kit



GSI Caspase-2 Colorimetric Assay Kit DataSheet

Assay Procedure

1. Treat cells by desired method. We recommend conducting four reactions (1) non-induced cells, (2) non-induced cells + inhibitor, (3) induced cells, and (4) induced cells + inhibitor as well as two control reactions (1) caspase 2 positive control, and (2) caspase 2 positive control + inhibitor.
2. Count cells and pellet 1×10^6 cells per test in 1.5 ml tube.
3. Resuspend cells in 100 μ l of chilled 1 x Assay Buffer and incubate cells on ice for 10 minutes.
4. Centrifuge at 10,000 x g for 1 min.
5. Transfer supernatant (cytosolic extract) to each well of 96-well plate (or a fresh tube if a spectrophotometer is later used) and put on ice for immediate assay (or store at -80°C for future use).
6. Immediately before use, prepare sufficient working reagents per assay: 10 μ l of Substrate, 5 μ l of 200 mM DTT and 85 μ l of 1 x Assay Buffer.
7. Add 100 μ l of the Working Reagents prepared as above to each well of 96-well plate.
8. Seal the plate with plate sealer. Incubate at 37°C for 1~2 hours and protect from light.
9. Read the plate at 405 nm (wavelength can range from 400 to 410 nm) in a microplate reader (or a spectrophotometer if microtubes are used).
10. Fold-increase in caspase-2 activities can be determined by comparing the OD_{405} with that of the uninduced control or other desired controls.

Precaution and Technical Notes

1. It is critical to follow the procedure step by step otherwise appropriate color development may not occur as expected.
2. A negative control without addition of substrate is recommended and the OD_{405} subtraction should be applied to each test to correct the background signal.
3. This kit should not be used beyond the expiration date on the label.
4. Use a fresh reagent reservoir and pipette tips for each step.
5. It is recommended that all controls and samples be assayed in duplicate.
6. Avoid microbial contamination of reagents and buffers. This may interfere with the sensitivity of the assay.